



Accessible. Durable.

AMSCO® Warming Cabinets

 STERIS

Options to Fit Your Needs

WARM, SECURE & ORGANIZED

AMSCO® Warming Cabinets are designed to store and warm sterile intravenous (IV) and surgical irrigation solutions, linens, and blankets to programmed temperature settings for patient applications. Cabinets are available in:

- Standard or premium featured configurations
- Single- or dual-compartment models
- Glass or solid door versions
 - 18" (457 mm) or 24" (610 mm) depths
 - Mobile or fixed-based models
 - Freestanding or recessed installations



Optional Mobile cart available for Single compartment model

Standard Configuration

PRECISE, CONVENIENT CONTROL

Our digital control technology offers:

- Pre-set temperature range for IV solutions, irrigation solutions, or blankets
- Incremental temperature adjustment capability
- Display shows either set-point or actual compartment temperatures in Fahrenheit or Celsius
- Visual indicator provides door ajar warning



Each compartment features an independent temperature control with pre-set buttons for IV temperature or IRR/blanket temperature.



Displays on dual compartment (left) and single compartment (above) can show temperature in F or C.

Premium Model Warming Cabinets

ADVANCED CARE WITH COMPLIANCE AND CONVENIENCE



Electronic data recording provides reports that include date and time, temperature setting, actual temperature, door status, over-temp status and under-temp status. Stores data up to 12 months.



On dual compartment models, the electric door lock feature will unlock both upper and lower compartments simultaneously. Key lock assures access in the event of a power failure.



Digital Data Recording and Keyless Locking

- Easy-access USB connection
- Automatic temperature monitoring, helps meet AORN guidelines
- Digital security to prevent inventory loss

Dual-pane tempered glass doors provide instant visual confirmation of supplies.

Mid-Size Single Compartment Warming Cabinet

PREMIUM CONTROLS

Digitally controlled warming in a size that fits your needs.

- Counter height unit can be freestanding, installed under a counter, or recessed
- Control panel includes electronic data recording and keypad door locking
- Optional mobile base (with bumpers) provides more flexibility
- 7.2 cubic feet of warm storage



Sturdy locking casters for mobility and cleaning

Two slide-out baskets for fluid storage are included with mid-size cabinet



Adjustable blanket storage shelf, also included, can hold blankets or bottles

Built-in full perimeter rubber bumpers protect unit (and facility walls)

Integrated Warming Unit

A SEAMLESS COMBINATION

AMSCO® Warming Cabinets combine with OR Storage Consoles for a sleek, integrated and organized storage space.



Pass-Thru Warming Cabinet provides the convenience of access inside and outside of the operating room, and gives staff constant inventory flow and direct access to supplies

Technical Specifications

Available Models	Single Compartment		Dual Compartment		OR Storage Warming Cabinet		Mid-Size Single Compartment	Pass-Thru Dual Warming Cabinet
	Standard	Premium	Standard	Premium	Standard	Premium	Premium	
Door Options	Steel or dual-pane tempered glass front				Steel front		Dual-pane tempered glass front	Steel front
Installation	Freestanding or recessed (Mobile cart available*)		Freestanding, recessed or mobile base (24" model)		Freestanding or recessed		Freestanding, recessed or mobile base	Full-recessed or partially-recessed
Outside Dimensions (H x W x D) 18" Depth	24¼ x 30 x 20½" (616 x 762 x 521 mm)		74¾ x 30 x 20½" (1899 x 762 x 521 mm)		60 x 30 x 18" (1524 x 762 x 457 mm)		n/a	60 X 24 x 18" (1524 x 610 x 457mm)
Outside Dimensions (H x W x D) 24" Depth	24¼ x 30 x 26½" (616 x 762 x 673 mm)		74¾ x 30 x 26½" (1899 x 762 x 673 mm)		n/a		36 3/8 x 30 x 27½" (924 x 762 x 689 mm)	n/a
Inside Dimensions (H x W x D) 18" Depth	13½ x 24 x 16⅞" (343 x 610 x 429 mm)		Upper Chamber 13½ x 24 x 16⅞" (343 x 610 x 429 mm)		n/a		n/a	Upper chamber 12¼ x 19 x 15¼" (311 x 483 x 387mm)
			Lower Chamber 36½ x 24 x 16⅞" (927 x 610 x 429 mm)		n/a		n/a	Lower chamber 24¾ x 19 x 15¼" (619 x 483 x 387mm)
Inside Dimensions (H x W x D) 24" Depth	13½ x 24 x 22⅞" (343 x 610 x 581mm)		Upper Chamber 13½ x 24 x 22⅞" (343 x 610 x 581 mm)		n/a		22¾ x 24 x 22⅞" (578 x 610 x 581 mm)	n/a
			Lower Chamber 36½ x 24 x 22⅞" (927 x 610 x 581 mm)		n/a		n/a	n/a
IV fluid- 510(k)	Yes							
Temperature range	90°F (32°C) to 160°F (71°C)							
Temperature Lockout	Yes							
Fahrenheit and Celsius	Yes							
Actual and set point temperature display	Yes							
Alarm	Audible and visual							
Alarm conditions	Over Temperature, Under Temperature and Door Ajar							
Door locking mechanism	Mechanical	Electronic with mechanical backup	Mechanical	Electronic with mechanical backup	Mechanical	Electronic with mechanical backup	Electronic with mechanical backup	Mechanical
Electronic Data Recording	No	Included	No	Included	No	Included	Included	No

*Mobile Cart (H x W x D) 37⁵/₈ x 36¹⁵/₁₆ x 28⁹/₁₆"

For more information, contact your **STERIS** representative or visit www.steris.com.